

Prelude™


Self-Etch, Total Etch & Dual-Cure Adhesive System

- ♦ **Fast & reliable** one 10 second application per bottle.
- ♦ Use of Prelude Link and self/dual-cure composites converts Prelude to dual-cure adhesive.
- ♦ **High bond strengths to uncut enamel!**
- ♦ Nano-filled, ethanol based formula is highly stable.
- ♦ Ultra thin film thickness, less than 5 micron.
- ♦ High 5 minute & 24 hour (34.7MPa) bond strengths.




Prelude, by Danville offers the latest in Self-Etch technology. One of the fastest and easiest to use, Prelude also offers high bond strengths to uncut enamel and self/dual-cure composites. Prelude Adhesive, used without Primer, also serves as a Total-Etch material. Nano-filler, ethanol-based formulation is highly stable and does not settle out.

Shear Bond Strength on Dentin at 5min. in self-etch technique



OptiBond Solo Plus, Clearfil SE, Unifill Bond, Xeno III are NOT registered trademarks of Danville Materials.

Fastest Application!



www.danvillematerials.com

3420 Fostoria Way, Suite A-200 San Ramon, CA 94583 800.827.7940 925.973.0710 Fax 925.973.0764



Innovative Design. Reliable Performance

92190 Rev B