



## Accolade SRO™


### Super Radiopaque Flowable Composite

- Superior radiopacity
- Superb esthetics
- Excellent polishability
- 0.7 micron fill
- Natural Fluorescence

With a radiopacity of 320% of aluminum, (similar to titanium) radiographs were never easier to read! Recently confirmed by a leading clinical reviewer, Accolade SRO is among the most radiopaque flowable on the market, yet is value priced. As a 0.7 micron hybrid microfill, Accolade SRO also offers exceptional strength, esthetics and the polishability of StarFlow and Accolade.

Shade	1gm	5gm
A-2	91415	91387
A-5	91416	91388
Extra Light	91417	91389

#### Radiopacity of popular flowable resins\*



\* Results by **Leading Clinical Reviewer**, call for details.

Esthet-X Flow, Filtek Supreme Plus Flow, Artiste, and Tetric Flow are NOT registered trademarks of Danville Materials.

[www.danvillematerials.com](http://www.danvillematerials.com)

3420 Fostoria Way, Suite A-200 San Ramon, CA 94583 800.827.7940 925.973.0710 Fax 925.973.0764



Innovative Design. Reliable Performance

92193 Rev B