

StarFlow™

Flowable Composite

StarFlow's superior flow characteristics and esthetics make ideal results easy. Sub-micron particle size gives superior polishability and resistance to staining. A leading clinical reviewer showed StarFlow to have the highest compressive strength of any flowable composite. True Vita* shade match, radiopaque and non-sticky. It offers fluoride slow release and natural fluorescence. (For higher viscosity, we offer Accolade)

StarFlow uses:


- Conservative Restorations
- Low wear applications
- Cavity Liner
- Refining TurboTemp 2. Ideal for use with all provisional composites

Intro Kit
 (3) 5gm syringes
 (A-2, A-3.5, & Incisal)
 (100) 20 gauge tips
 Ref # 85071



Trial Kit
 (4) 1gm syringes
 (A-1, A-2, A-3.5, & Incisal)
 (20) 20 gauge tips
 Ref # 89485

Compressive Strength of Flowable Resins



Individual Syringes					
Shade	1gm	5gm	Shade	1gm	5gm
A-00op	89465	89421	A-5	89471	85041
A-0	89466	89397	B-1	N/A	90471
A-1	89467	85051	B-2	89472	85055
A-2	89468	85052	C-2	89473	85056
A-3	89469	85053	C-4	89474	85057
A-3.5	89470	85054	Dentin	89476	85058
A-4	90749	90750	Incisal	89475	85059

Aria™

Microfill Flowable Composite

Aria's microfill formulation offers polishability and flexibility no hybrid can match! Ideal for facials and Class V. High translucency gives Aria "chameleon" qualities, necessitating only a limited number of shades. Aria seems to disappear into the surface making it an ideal companion to Accolade and StarFlow when less opacity is desired.

Shade	1gm	5gm
A-1	89356	89360
A-2	89357	89361
A-3.5	89358	89362
Incisal	89359	89363

Starter Kit
 1gm of each shade &
 (20) 18 gauge tips
 Ref # 89355



* Not a trademark of Danville Materials

www.danvillematerials.com

3420 Fostoria Way, Suite A-200 San Ramon, CA 94583 800.827.7940 925.973.0710 Fax 925.973.0764



Innovative Design. Reliable Performance.

92194 Rev A