

Intraoral Sandblasting

Sandblasting is fast and easy with most surface roughening procedures requiring only a few seconds. All models feature carbide tips and two abrasive jars for rapid abrasive change. (ask Craig) All MicroEtchers include a 1 year warranty.

Each model includes:

- Stainless steel handpiece
- .048 carbide nozzle assembly
- 9' air line
- Instruction manual



MicroEtcher™ IIA

Ref # 22005-01

- The MicroEtcher IIA has a permanently affixed jar lid and spill resistant abrasive filter.
- Fully Autoclavable 360° rotatable, slim profile nozzle offers better intraoral usability.

MicroEtcher™ II

Ref # 22000-03



- Only the nozzle is autoclavable.
- Replaceable carbide tip.

Optional Nozzles Sandblaster Nozzles

Interchangeable nozzles for the MicroEtcher II and IIA. New smaller .032" tip is ideal for intraoral applications (pit & fissure preps or porcelain repair) as it reduces the amount of dust and provides for a more focused pattern of abrasive.

II Nozzle
.048 60°
Ref # 201802-00

Pit & Fissure Nozzle
.032 60°
Ref # 186113

**Slim Profile
IIA Nozzle**
.048 60°
Ref # 201865-00



MicroEtcher Uses:

- Etching all metals, composites & amalgam for maximum bond strength
- Aiding in intraoral porcelain repair
- Prepping for pit & fissure restoration
- Cleaning dark grooves to reveal decay
- Prepping incipient lesions and minor Class I & V
- Debriding bur preps for best esthetics & bonds

www.danvillematerials.com

3420 Fostoria Way, Suite A-200 San Ramon, CA 94583 800.827.7940 925.973.0710 Fax 925.973.0764



Innovative Design. Reliable Performance

92196 Rev A