

SpinBright™

DIAMOND POLISHING BRUSHES

- ◆ POLISH COMPOSITE AND PORCELAIN
- ◆ DIAMOND IMPREGNATED BRISTLES OFFER SUPERIOR POLISH AND LIFE
- ◆ AUTOCLAVABLE FOR UP TO 50 USES

Diamond impregnated bristles polish composites and porcelain rapidly. Diamond offers superior results and durability compared to competitors', carbide brushes. Autoclavable for as many as 50 uses. Latch type, 2 shapes available. 10 brushes per pack.



SpinBright Intro Kit 10 brushes (5 universal, 5 occlusal)	90917
SpinBright Universal 10 universal brushes	90919
SpinBright Occlusal 10 occlusal brushes	90918



- ◆ CREATE GREATER ACCESS TO THE SULCUS WHEN TAKING IMPRESSIONS.
- ◆ REDUCE GINGIVA FOR RESTORING DEEP CLASS V CAVITIES.
- ◆ TRIM TISSUE AFTER RECOVERY FROM IMPLANT SURGERY.

Trim Eze **91200**
One Bur

TrimEze™

HIGH SPEED TISSUE TRIMMER

TrimEze produces similar results to electro surge or laser without the cost! Narrow tip for good access. Minimally invasive that promotes healing and regrowth. Narrow tip is convenient even in areas of difficult access. Use in your high speed handpiece without water spray.

TrimEze can be autoclaved for reuse. Expected life is up to several hundred uses.

www.danvillematerials.com

3420 Fostoria Way, Suite A-200 San Ramon, CA 94583 800.827.7940 925.973.0710 Fax 925.973.0764



Innovative Design. Reliable Performance

93142 Rev A