

Retract™

GINGIVAL RETRACTORS



TWO SIZES AVAILABLE

NON-SLIP GRIP

AUTOCCLAVABLE

CONTOURED BLADES



Retract, Standard
Retract, Large

89445
89448

Retract is ideal for gingival retraction during Class V preparation and for veneer placement. The blades are sharpened to allow composite shaping. Retract offers left and right hand blades with blunted points to prevent gingival injury. Two blade widths are offered, standard and large. The large size is also ideal for retracting and shielding during Class V air abrasion.

Retract is made of surgical stainless and is fully autoclavable.

www.danvillematerials.com

3420 Fostoria Way, Suite A-200 San Ramon, CA 94583 800.827.7940 925.973.0710 Fax 925.973.0764



Innovative Design. Reliable Performance

93146 Rev A