


Danville's *StarFill 2B*™



The First Dual-Cure, Automix, Directly Injectable Composite!

The same great product, now with a fast & easy delivery system!
Stop wasting time hand mixing Base & Catalyst. With StarFill 2B's new Dual Automix Syringe, direct placement is easy.



Features:

- Dentin Shade
- Dual-Cure
- Fluoride Slow Release
- High Strength & Polishability
- More Economically Priced!
- Bonds to Both Light-Cured & Auto-Cured Composites

Ideal Luting System

StarFill 2B is ideal for luting posts & crowns where light curing is not possible. StarFill 2B adapts easily, making it ideal as a cavity liner and assuring full polymerization even in undercuts. Dual-cure chemistry allows for either rapid light-cure or self-cure in 3 min. Both light-cure and self-cure composites can be placed directly over StarFill 2B.

Directed Shrinkage

This results in reduced post-op sensitivity and enhanced bond strength. When used as a cavity liner StarFill 2B tends to shrink towards that tooth.

StarFill 2B
(10gm syringe & 6 automix tips)

90150

Automix Tips (12)
Automix Tips (48)

90366
90368

www.danvillematerials.com

3420 Fostoria Way, Suite A-200 San Ramon, CA 94583 800.827.7940 925.973.0710 Fax 925.973.0764



Innovative Design. Reliable Performance

93290 Rev A