


A major breakthrough in air abrasion systems!




PREPAIR™

Pit & fissure & etching handpiece

The PrepAir is a handheld air abrasion unit with unmatched performance. The finger button control gives instantaneous on/off response and eliminates the need for a foot pedal. Abrasive flow is even and set to a minimum, to greatly reduce over-spray. Interchangeable nozzles rotate 360° assuring proper reach from any angle. Hook up is easy as the PrepAir will connect to a MicroEtcher™ disconnect, high-speed outlet or numerous other air connections.



"Fissure" Bur



PrepAir

PREPAIR

APPLICATIONS AND BENEFITS

- Clean darkened grooves to reveal decay
- Clean for reliable Diagnodent reading
- Most conservative & precise
- Small cavity prep
- Pit & fissure prep
- Etch orthodontic appliances
- Etch restorations for maximum bond
- Etch restorations for maximum bond
- Etching for lingual retainers

www.danvillematerials.com

3420 Fostoria Way, Suite A-200 San Ramon, CA 94583 800.827.7940 925.973.0710 Fax 925.973.0764



Innovative Design. Reliable Performance

93291 Rev A