

# Ice Light™

## Light Transmitting Post



**ICE**LIGHT is the next generation of aesthetic posts from Danville! **ICE**LIGHT features 70% fill of Zirconia enriched, parallel glass fibers for exceptional strength. Post design is parallel with a tapered end. Its flexibility is close to dentin so stress is distributed, minimizing the risk of cervical fractures.

**ICE**LIGHT is light transmitting for use with light-cured as well as self/dual-cured composites and adhesives. Posts are also presilanated, making placement rapid. The length is easily trimmed eliminating the stocking of multiple lengths. One reamer, per size, gives the ideal post clearance

Non-light transmitting, **ICE**POST, are also available. Both are value priced.

and allows excess composite to escape. Further, there is no composite post with higher radiopacity!

**ICE**LIGHT kit feature 20 posts in sizes 1mm, 1.2, 1.4 & 1.6mm with reamers.



**ICEPOST Intro Kit 20/pk 90079-01**  
 Econo Kit 20/pk 90084-01  
 1.0 mm 10/pk Refill 90082-01  
 1.2 mm 10/pk Refill 90083-01  
 1.4 mm 10/pk Refill 90088-01  
 1.6 mm 10/pk Refill 90089-01

**ICE**LIGHT Intro Kit 20/pk **90080-01**  
 1.0 mm 10/pk Refill 90378-01  
 1.2 mm 10/pk Refill 90380-01  
 1.4 mm 10/pk Refill 90382-01  
 1.6 mm 10/pk Refill 90384-01

[www.danvillematerials.com](http://www.danvillematerials.com)

3420 Fostoria Way, Suite A-200 San Ramon, CA 94583 800.827.7940 925.973.0710 Fax 925.973.0764



Innovative Design. Reliable Performance

93294 Rev A