

# MicroPrime B™

## Universal Desensitizing Agent

- Quick & Easy Application
- Universal Compatability
- Antimicrobial Agent
- Will Not Burn Soft Tissue



Eliminate  
Post-op  
Sensitivity!  
*without soft tissue irritation.*

**B** MicroPrime is a superior desensitizing agent which can be placed under dental cements or temporary, provisional or final restorative materials to reduce post-operative sensitivity.

It can also be used for desensitization of amalgam restorations, either conventional or bonded. MicroPrime contains HEMA to aid bonding primers to penetrate etched dentin. The HEMA base in MicroPrime helps seal tubules, giving immediate relief from pain due to hypersensitivity.

**APPLICATION IS QUICK AND EASY.**

**MICROPRIME IS IDEAL FOR:**

- Abrasion, cervical erosions and preps
- Applications to vital crown preparations prior to luting with cements.
- At the “prep” appointment to desensitize during temporization.

Most dentin resin bonding materials will benefit from MicroPrime application. Its application reliably reduces post-op sensitivity by supporting the collagen framework for easier penetration of the adhesive, thus enhancing the dentin bond.

[www.danvillematerials.com](http://www.danvillematerials.com)

3420 Fostoria Way, Suite A-200 San Ramon, CA 94583 800.827.7940 925.973.0710 Fax 925.973.0764



Innovative Design. Reliable Performance

93296 Rev A