

## Microfill Flowable Composite

- Ideal for:
- Enamel Only Restorations
- All Class V
- Incisal Edge Repair
- As a Final Layer in Composite restorat



### Features:

- Fast & Easy.
- Highly filled (61% filled).
- Excellent access via the syringe.
- Microfiller Particle size = 0.04 microns.
- Available in selected vita shades.
- Superior translucency and handling.
- Particularly well suited for Air Abrasion.
- Uniform submicron particle size.

### Benefits:

- Injected directly into the preparation.
- Wear rate similar to natural tooth structure.
- Facilitates fast and easy placement.
- Excellent strength and esthetic.
- Allows for a variety of cases and situations.
- “Chameleon” effect & will not slump.
- Great flow characteristics.
- Excellent marginal adaptation.
- Extreme polishability &

**A**ria is a light-cured, microfill flowable composite, designed for use as the “final or enamel layer” in any composite restoration.

Aria is so lifelike that it resembles natural tooth structure, combining superior translucency with high strength, better esthetics and greater wear properties. It has the superior polishability of a microfill, and the ease of placement and the versatility of modern flowable composites. It offers high flow and superior handling characteristics. The viscosity provides for excellent marginal adaptation.

**89355 Aria Introductory Kit**  
4 - 1gm syringes  
(A-1, A-2, A-3.5, Incisal)  
20 - 18 gauge tips

**89356** Aria 1g A-1  
**89357** Aria 1g A-2  
**89358** Aria 1g A-3.5  
**89359** Aria 1g Incisal

**89360** Aria 5g A-1  
**89361** Aria 5g A-2  
**89362** Aria 5g A-3.5  
**89363** Aria 5g Incisal

Aria 1 and 5 gram refills come with 6 - 18 gauge tips.

