

Caries Finder™

Caries Disclosing Dye

- Prevent Recurrent Caries
- Retain Maximum Tooth Structure
- Available in Red or Green



Caries Finder helps retain maximum tooth structure and reduce possibility of leaving carious dentin .

By precisely staining only carious dentin or enamel Caries Finder dye ensures conservative caries removal. It reduces the possibility of inadvertently leaving carious dentin or of removing excessive dentin. Caries Finder will not stain the following materials:

- Normal healthy dentin.
- Inner carious dentin (uninfected) with solid collagen fibers.
- Normal healthy enamel.

Since Caries Finder only stains the denatured collagen fibers of infected dentin, healthy and remineralizable dentin is easily differentiated. It saves time by identifying incipient lesions, dentin caries, decay in undercuts and lesions with limited view or access. It is especially useful where pulp exposure is expected by indicating where excavation should stop.

Caries Finder Green 10ml 80004
Caries Finder Red 10ml 80010
Caries Finder Red 1gm 90764
Syringable

QUICK AND EASY APPLICATION

