

# ABC WEDGE™

ABSOLUTE • BEST • CONTACTS

## Sectional Matrix System

- Restores natural contour
- Ideal contacts
- Minimizes finishing
- Retentive Tines.
- Flat Tines.
- Contoured stainless steel Matrices.
- Occlusal offset.
- Stackable.

...its as easy as

# 1-2-3



ABC Wedge™ with Contact Matrix™ is the ultimate sectional matrix system, with three improvements over any other matrix system. The ABC Wedge makes restoring the wide open box or missing cusp easy when using sectional matrix or trimmed Tofflemeyer bands. It is all too frequent that excessive decay results in fractured cusps. Previously, the only alternative was a Tofflemeyer band with the obtrusive retainer, frequently creating difficulty in restoring a tight, biomimetic contact.

**ABCWedge Kit Ref # 93488**  
 (192 Wedges)  
 X-Small (24R, 24L)  
 Small (24R, 24L)  
 Medium (24R, 24L)  
 Large (24R, 24L)  
 Mega Rings (2) Gold PN 91304  
 Pliers (1) PN 91298

**Contact Rings feature:**

1. **Converging tines** to make the rings far more retentive than parallel tines.
2. **Inward and outward tines** for placement on either side of the cavity preparation.
3. **Occlusal offset** to allow stacking and also clearance over obstacles like rubber dam clamps.

**Refills (192 Wedges)**  
 X-Small Ref # 93486  
 Small Ref # 93487  
 Medium Ref # 93588  
 Large Ref # 93589



1 Failed Amalgam



2 Fractured Buccal Cusp



3 ABC Wedge Supports the matrix



4 Final Restoration

[www.danvillematerials.com](http://www.danvillematerials.com)

3420 Fostoria Way, Suite A-200 San Ramon, CA 94583 800.827.7940 925.973.0710 Fax 925.973.0764



Innovative Design. Reliable Performance

93619 Rev A