

# AlphaBite™

## Bite Registration Material

Alpha Bite is Danville's new bite registration material that is formulated to the specifications of Blu Mousse®<sup>1</sup>. It has the same 'standby' properties, yet is more value priced. It is also easy to carve without breaking. Alpha Bite is made from quality fillers and silicones. Work time is 45 sec, Set is 75 sec. Available in 50ml cartridges.

<sup>1</sup>Blu Mousse is a trademark of Parkell.



**AlphaBite**  
(2) 50ml + 6 tips  
Ref # 93352



**Noginate**  
(2) 50ml + 6 tips  
Ref # 93553

## Noginate™

### VPS Alginate Substitute

Noginate is a value priced VPS material that may be used as an alginate substitute. Ideal for producing orthodontic models, study models, for creating mouth guards and bleaching trays. Unlike alginate, Noginate is very stable, has a high tear strength and may be sterilized. Noginate is so stable that it can be archived instead of heavy stone models. Work time: 90 sec. set: 150 sec. Value priced, available in 50ml cartridges.

[www.danvillematerials.com](http://www.danvillematerials.com)

3420 Fostoria Way, Suite A-200 San Ramon, CA 94583 800.827.7940 925.973.0710 Fax 925.973.0764



Innovative Design. Reliable Performance.

93620 Rev A